

WARRANTY

STABILA[®] TYPE 187 K

Spirit Levels

STABILA spirit levels will maintain their stated accuracy for life. The vials on all STABILA levels will not fog, leak, or become inaccurate. If they do, you will be given a new level. **Bent frames** void warranty.

Laser and Tech (Electronic) Levels

We will repair, without charge, any defects due to faulty materials or workmanship for two years from the date of purchase.

Accessories: Tape Measures, Tripods, Cases, Grade Rods, Etc.

30 Day Warranty against manufacturer defects from the date of purchase.

Is The Level Accurate?

To check the "level" or "horizontal" vial's accuracy, place the level on any flat surface.

- + Make a mark on the flat surface at one end of the level.
- + Read the bubble & remember where it is.
- + Turn the level 180 degrees (end to end), placing the level right on the line drawn from the original end placement.
- + Read the level.
- + If the bubble returns to the same place, the level is accurate. If it does not, it is not accurate.
- + Do it again, regardless of the reading.

To check the "plumb" or "vertical" vials, place the level against any door frame or wall.

- + Make a mark on the wall beside the top vial.
- + Read the bubble & remember where it is.
- + Flip the level side for side (not end for end!).
- + Read the level.
- + If the bubble returns to the same place, the vial is accurate. If it does not, it is not accurate.
- + Do it again, regardless of the reading.

To check the second "plumb" vial, turn the level end for end.

- + Repeat all of the steps as done for the first "plumb" vial.

Laser Levels should always be returned to Stabila for warranty evaluation.

All Stabila returns require a Returned Goods Authorization (RGA). Please contact our Customer Service Department to have an RGA issued.

You can also visit your nearest Stabila Distributor to process your warranty return.